

Color shades

Kuraray Dental Color Code

US	Universal Standard	DG	Dark Grey	YBO	Yellow Brown Opaque
UY	Universal Yellow	LO	Light Opaque	DGO	Dark Grey Opaque
UL	Universal Light	YO	Yellow Opaque	L	Light
UO	Universal Opaque	HO	Hollywood Opaque (Bleach White)	U	Universal
XL	Extra Light	CL	Clear Transparent	Y	Yellow
DY	Dark Yellow	DYO	Dark Yellow Opaque	G	Grey
YB	Yellow Brown			B	Brown

Conversion of the Vita® Color Code

Vita	A1	A2	A3	A3.5	A4	B1	B2	B2.5	B3	B4	C2	C3
*KURARAY	UL	US	DY, (B)	G	–	L, XL	U	UY	Y, YO	YB	DGO	DG

Vita	OA1	OA2	OA3	OB4	Enamel	HO
*KURARAY	LO	UO	DYO	YOB	–	–

*For the following products only:
CLEARFIL™ PHOTO POSTERIOR, PHOTO CLEARFIL™ BRIGHT,
CLEARFIL™ ST OPAQUER, CHROMA ZONE™ COLOR STAIN.

Conversion of the Kuraray Dental Color Code to the Vita® Color Code

Kuraray Color Code	Color description	Shade	Vita Color Code
Enamel	Enamel		Enamel
HO	Hollywood Opaque (Bleach White)	Light	} A1
UL	Universal Light	Dark	
LO	Light Opaque	Opaque	
US	Universal Standard		} A2
UO	Universal Opaque	Opaque	
DY	Dark Yellow	Light	} A3
B	Brown	Dark	
DYO	Dark Yellow Opaque	Opaque	
G	Grey		A3,5
XL	Extra Light	Light	} B1
L	Light	Dark	
U	Universal		B2
UY	Universal Yellow		B2,5
Y	Yellow		} B3
YO	Yellow Opaque	Opaque	
DYO	Dark Yellow Opaque	Opaque	
YB	Yellow Brown		} B4
YBO	Yellow Brown Opaque	Opaque	
DGO	Dark Grey Opaque	Opaque	C2
DG	Dark Grey		C3